

DERMA
HOME



UNI&CORE

DERMA HOME

DERMA HOME Guidebook

ALLTHERMA™

Built on the RF principle used in dermatology clinics, 1,000,000 Hz radio-frequency energy reverses direction a million times per second to reach deep into the skin. 10-level RF care focused on sagging along the facial contours, jawline, and the eye and mouth areas.




How-to Video QR

Recommended level 2-3

Recommended time 10 min

Scan the QR code with your camera to watch the how-to video.


Areas Face / Jaw / Shoulders



From deep within,
smooth and resilient —
powerful tightening

ULTRASONIC

Using the ultrasound principle of clinical skin care, 800 kHz–1 MHz energy vibrates up to one million times per second and travels deep into the skin. It stimulates collagen and elastin production to restore firmness and volume to sagging cheeks, nasolabial folds and the jawline.



How-to Video QR

Recommended level 2–3

Recommended time 10 min

Scan the QR code with your camera to watch the how-to video.

Areas Cheeks / Nasolabial folds / Deep wrinkles

Sagging cheeks,
deepening smile lines —
Ultrasonic restores
firmness and volume

IONTOPHORESIS

The same 2-way iontophoresis principle used in clinics. It cleanses impurities from deep within the pores, then delivers ionized active ingredients deep into the skin for absorption far beyond simple application.



How-to Video QR

Recommended level 2-3

Recommended time 10 min

Scan the QR code with your camera to watch the how-to video.

Areas Face / Neck / Shoulders

Advanced technology,
maximized absorption —
the **intuitive solution**
for skin concerns




Home Screen

Touch the screen
to go to the next step.



Language

Select your language
(ENGLISH).



Select Function

Touch the function
you want to use.



Apply Care Product

Apply the UNI&CORE product
you wish to use to your face.




Function Detail

Set the time and intensity
with the buttons, then start.

Firming Lift Routine **RF Radio Frequency**

An intensive solution to lift sagging skin back up


- ① Apply **UNI&CORE Deep Repair Nutrition Cream** over the entire face.
- ② In **RF Radio Frequency**, set level 2–3 and 10 minutes, then press Start.
- ③ Glide the head from under the chin along the jawline toward the ears, lifting upward.
- ④ Repeat on the cheeks and forehead, always moving upward and outward.
- ⑤ Finish by gently patting any remaining cream into the skin.



Radiance Anti-aging Routine **Iontophoresis**

A clear-skin solution that maximizes the power of Vitamin C


- ① Apply **UNI&CORE Vita C 13.5 Whitening Ampoule** over the entire face.
- ② Hold the **Deep Clean (A)** handle and select **Iontophoresis**.
- ③ Set level 2–3 and 10 minutes, then press Start.
- ④ Sweep the **Cleansing (N)** head gently over the skin to draw out impurities.
- ⑤ Switch to the **Absorption (C)** head and sweep the same areas again to drive in the actives.
- ⑥ Finish with **UNI&CORE Deep Repair Nutrition Cream**.



Face All-care Routine **Ultrasonic + RF + Iontophoresis**

Where care begins — a total solution for troubled skin

- ① Apply **UNI&CORE Soothing & Moist Care Boosting Gel** generously over the face.
- ② In **Ultrasonic**, set level 2–3 and 10 minutes, then press Start.
- ③ Wipe off the gel and cleanse the device head.
- ④ Follow steps 1–5 of the **Radiance Anti-aging Routine**.
- ⑤ Finish with the **Firming Lift Routine** using Deep Repair Nutrition Cream.
- ⑥ Alternate the Firming and Radiance routines for 1–2 weeks for best results.



Can I use the device every day?

Yes. Each function is designed for daily use within the recommended time and intensity. If your skin feels sensitive, rest for a day and resume at a lower level.

What should I apply before use?

Always use a UNI&CORE gel, ampoule or cream matched to the function. Using the device on dry skin reduces conduction and may cause discomfort.

Is it OK to raise the intensity?

Start at level 2–3 and increase gradually only once your skin is comfortable. Higher levels do not shorten the time needed for results.

Can I use other brands' products?

Water-based products are generally compatible, but oil-based or exfoliating products should be avoided. Results are optimized with the paired UNI&CORE line.

Precautions

A. Do not use if you

- Have a pacemaker, implanted electronic device, or metal implants in the treatment area
- Are pregnant or breastfeeding, or have epilepsy or a heart condition
- Have skin disease, open wounds, inflammation, eczema or severe acne on the area
- Have had filler, Botox or laser treatment in the area within the last 2 weeks
- Have a malignant tumor, diabetes, thrombosis or sensory impairment of the skin
- Are under 18 years of age or unable to express discomfort
- Have a fever, infectious disease, or are under the influence of alcohol

B. Product specifications

- Product name: Home Beauty Device
- Model: Derma Home 10 Pro
- Rated voltage: AC 110–250 V / 50–60 Hz
- Size / weight: 29.5 (W) × 29.5 (D) × 115.5 (H) cm / 19 kg
- Manufacturer / Country of origin: UNI&CORE Co., Ltd. / Republic of Korea

C. Handling precautions

- Read this guide carefully before use and keep it for future reference.
- Use only the supplied power cord and connect to a grounded outlet.
- Keep the device away from water, steam and direct sunlight; never use in the bathroom.
- Do not use on the eyelids, lips, thyroid area or over the heart.
- Do not exceed the recommended time and intensity for each function.
- Stop immediately if you feel pain, burning or redness that does not subside.
- Clean the heads with a soft, dry cloth after each use. Do not use alcohol or abrasive cleaners.
- Do not disassemble, modify or repair the device yourself.
- Keep out of reach of children; this device is not a toy.
- Use only UNI&CORE care products designed for this device.
- Unplug the device when not in use and before cleaning.
- Contact customer service if the device is dropped, damaged or operates abnormally.

UNI&CORE

Warranty Certificate

This product has passed strict quality control and is warranted in accordance with the consumer protection regulations of the Republic of Korea.

Free repair is provided for manufacturing defects occurring under normal use within the warranty period.

* If a defect occurs during normal use within the warranty period, it will be repaired free of charge under these warranty terms.

The warranty period is 1 year from the date of purchase.

* Damage caused by misuse, negligence, unauthorized repair, natural disaster or normal wear is excluded from free service.

Repairs after the warranty period or for excluded causes are provided at cost.

* For repair requests or inquiries, contact the UNI&CORE A/S Center at +82-2-571-0113 or your place of purchase.

Product		Model	
Installation address			
Dealer		Purchase date	
Customer	Name	Phone	
	Address		

Best paired with your device

SKINCARE

DERMA HOME



Scan the QR code with your camera
to see product details.



Deep Repair Activating Nutrition Cream 30 ml

A real skin-firming cream with powerful actives in a highly absorbable formula. Fortified with repairing ceramides and peptides, it restores elasticity to tired skin and reinforces the moisture barrier, leaving skin smooth and resilient.

Vita C 13.5 Super Whitening Ampoule 30 g

A 7-day brightening ampoule with 13.5% live, pure Vitamin C. Its self-activating formula targets dullness and dark spots the moment it meets the skin, for a visibly clearer, more even tone.

* Functional brightening cosmetic. ** Non-irritation tested.



Vita C Super Brightening Mask 5 Sheets

A radiance mask with 2% live, pure Vitamin C. Highly concentrated essence soaks into a fine sheet that clings closely to the skin, brightening and refreshing in a single 10-minute application.

* Functional brightening cosmetic.

Soothing & Moist Care Boosting Gel 250 ml

Skin-boosting care in a moist gel texture with patented botanical ingredients. Calms, hydrates and improves conduction for Ultrasonic care, so energy reaches the skin evenly without drying or irritation.

Gentle texture + light application + lasting hydration